


# RaySafe i3

– with on-screen display


INSTRUCTIONS FOR USE (EN) .....	3
BRUGSANVISNING (DA).....	4
GEBRAUCHSANWEISUNG (DE) .....	5
INSTRUCCIONES DE USO (ES).....	6
KÄYTTÖOHJEET (FI) .....	7
MODE D'EMPLOI (FR).....	8
ISTRUZIONI PER L'USO (IT) .....	9
使用説明書 (JA).....	10
GEBRUIKSAANWIJZING (NL) .....	11
BRUKSANVISNING (NO) .....	12
INSTRUÇÕES DE USO (PT).....	13
BRUKSANVISNING (SV) .....	14
使用说明 (ZH).....	15
APPENDIX.....	16

2.

3.

4.

# INSTRUCTIONS FOR USE (EN)

## Intended use

RaySafe i3 is designed for personal use on hospital staff. It measures scattered radiation in examination rooms using a dosimeter, and shows a graphical visualization of dose rate in real time on a display.

## Usage

The RaySafe i3 system needs to be installed and initiated before use. This can be done by a service technician, using the instructions in the printed *Installation and Service Manual*, delivered with the system.

Use the dosimeter on the outside of your X-ray protection garment. Front forward, clip towards your body.

The Real-time Dosimeter (**A** in Figure "1." on page 2) measures scattered radiation and stores measurements. When a dosimeter is exposed to radiation while connected to a Real-time Hub (**B** in Figure "1." on page 2), measurements are transferred from dosimeter to hub via radio. A display shows dose rate bars in real time in green, yellow or red, from low to high dose rate (**C** in Figure "1." on page 2).

To view dose history stored in the dosimeter, connect the dosimeter to a computer running Dose Viewer (**D** in Figure "1." on page 2). See the *Dose Viewer manual*, found in the help menu in Dose Viewer. Optional analysis software: Dose Manager (**E** in Figure "1." on page 2).

**CAUTION!** *The Real-time Dosimeter measures scattered radiation and should neither be positioned in the primary X-ray field nor be used for diagnostic measurements.*

**WARNING!** *RaySafe i3 does not replace any X-ray protection equipment.*

**WARNING!** *RaySafe i3 is not sterile.*

**WARNING!** *RaySafe i3 is not a replacement for legal dosimeters.*

## Maintenance

The Real-time Dosimeter has a replaceable battery. When a red light is flashing on the side of the dosimeter, the battery needs to be replaced. A service technician can do this, using the instructions in the *Installation and Service Manual*.

For cleaning of the dosimeter, use a damp cloth and mild detergent.

## Specifications

### Physical

	Real-time Dosimeter	Real-time Hub	Video Unit
Weight	34 g (1.2 oz)	180 g (6.3 oz)	355 g (12.5 oz)
Dimensions	40 × 58 × 17 mm (1.6 × 2.3 × 0.7 in)	184 × 134 × 35 mm (7.2 × 5.3 × 1.4 in)	112 × 84 × 34 mm (4.4 × 3.3 × 1.3 in)
Power	Replaceable battery: CR2450 Li/MnO <sub>2</sub>	Power over Ethernet: GT-963003648-T3-PP	External supply: EA1024PR
Video output			HDMI 1.4b

### Radiology

Operational quantity	H <sub>p</sub> (10)
Energy dependence	< 25 % (N-series, 40 – 150 kV)
Temperature dependence	< 5 % (18 – 26 °C) < 25 % (15 – 18 °C, 26 – 35 °C)
Response time	< 1 s (> 100 µSv/h), < 5 s (< 100 µSv/h)
Integration interval	1 s
Detection limit	< 30 µSv/h
Dose rate uncertainty (continuous radiation)	10 % or 10 µSv/h (40 µSv/h – 150 mSv/h), 20 % (150 – 300 mSv/h), 40 % (300 – 500 mSv/h)
The measured dose rate is monotonically increasing up to 1 Sv/h.	
Dose reproducibility	10 % or 1 µSv
Maximum lifetime dose	10 Sv
Typical angular response (relative to 0°)	Vertical: Figure "2." on page 2 Horizontal: Figure "3." on page 2
Typical energy response	Figure "4." on page 2
Real-time Hub dose unit	Sv or rem

## Radio

	Real-time Dosimeter	Real-time Hub
Frequency	866.0 MHz, 868.3 MHz, 918.3 MHz or 927.9 MHz, depending on purchased configuration (see label on product)	
Transmitter power output	3.2 mW	1 mW
Communication range	Typically 5 – 15 m, depending on the local environment	

### Climatic conditions

Temperature	+15 – +35 °C (for operation) –25 – +70 °C (for storage)
Temperature rate of change	< 0.5 °C/minute
Relative humidity	20 – 80 % (for operation) 5 – 95 % (for storage)
IP code	Real-time Dosimeter: IP20 Real-time Hub: IP20

## For safe usage

Always power off the dosimeter before air transport, to make sure that the radio communication is completely turned off.

## Security and privacy

Dosimeter serial number as well as name and settings entered via Dose Viewer are sent unencrypted via radio to the Real-time Hub.

## Symbols on labels

Symbol	Meaning
	Manufacturer
	Date of manufacture
	Article number
	Serial number
	Consult the instructions for use.
	Conforms to European Union directives.
	UK Conformity Assessed.
	Electrical safety compliance marking, together with listing number ("Appendix" on page 16)
	Radio certification number for US market ("Appendix" on page 16)
IC	Radio certification number for Canadian market ("Appendix" on page 16)
	Radio certification symbol and number for Brazilian market ("Appendix" on page 16)
	Conforms to relevant Australian Safety, radio and EMC standards..
R-NZ	Conforms to relevant New Zealand radio standards.
	Radio certification symbol and number for Japanese market
	In the European Union the WEEE directive and label indicates that when discarding this product, it must be sent to separate collection facilities for recovery and recycling. For advice and information, contact RaySafe.
	None of the China RoHS restricted substances are present above permitted levels.

## Learn more

- Real-time Dosimeter quick guide under the lid of the dosimeter packaging
- Dose Viewer manual via the Help button in Dose Viewer
- Installation and Service manual delivered with the system

## Service and support contacts

For information about service and support, visit [www.raysafe.com](http://www.raysafe.com).

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# APPENDIX

## MET

MET Laboratories, Inc. Certification covers UL 60950-1/CSA C22.2 No. 60950-1. MET Laboratories, Inc. has not evaluated this device for reliability or efficacy of its intended functions.

## FCC

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

*NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.*

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This portable device with its antenna is excluded from SAR testing since the output power is below the SAR test exclusion threshold. The antenna used for this device must not be co-located or operating in conjunction with any other antenna or transmitter.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

## IC/ISED

This Device complies with Industry Canada License-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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1256-14-5254

## Declaration of Conformity

Hereby, Unfors RaySafe AB declares that the radio equipment type RaySafe i3 is in compliance with Directive 2014/53/EU and UK Radio Equipment Regulations 2017.

The full text of the EU Declaration of Conformity and the UK Declaration of Conformity can be downloaded from the RaySafe i3 product page on [www.raysafe.com](http://www.raysafe.com)

## EU Import

EU importer of GlobTek power over ethernet GT-96300-3648-T3-PP:

Unfors RaySafe AB  
Uggledalsv. 29  
427 40 Billdal, Sweden

EU importer of EDAC power supply EA1024PR:

System Technology Sweden AB  
Linnégatan 14  
582 25 Linköping, Sweden

## Sustainability

For declaration of regulated substances, see <https://www.raysafe.com/about-us/sustainability>