
RaySafe i3



– with stand-alone display

INSTRUCTIONS FOR USE (EN)	3
BRUGSANVISNING (DA).....	4
GEBRAUCHSANWEISUNG (DE)	5
INSTRUCCIONES DE USO (ES)	6
KÄYTTÖOHJEET (FI)	7
MODE D'EMPLOI (FR)	8
ISTRUZIONI PER L'USO (IT)	9
使用説明書 (JA).....	10
GEbruiksaanwijzing (NL)	11
INSTRUÇÕES DE USO (PT)	12
BRUKSANVISNING (SV)	13
使用说明 (ZH).....	14
APPENDIX.....	15
EU DECLARATION OF CONFORMITY	16


1.


2.

3.

4.

5.

INSTRUCTIONS FOR USE (EN)

Intended use

RaySafe i3 is designed for personal use on hospital staff. It measures scattered radiation in examination rooms using a dosimeter, and shows a graphical visualization of dose rate in real time on a display. (Figure “1.” on page 2.)

Usage

The RaySafe i3 system needs to be installed and initiated before use. This can be done by a service technician, using the instructions in the printed *Installation and Service Manual*, delivered with the system.

Use the dosimeter on the outside of your X-ray protection garment. Front forward, clip towards your body.

The Real-time Dosimeter (**A** in Figure “2.” on page 2) measures scattered radiation and stores measurements. When a dosimeter is exposed to radiation while in range of a Real-time Display (**B** in Figure “2.” on page 2), measurements are transferred from dosimeter to display via radio. The display shows dose rate bars in real time in green, yellow or red, from low to high dose rate. You can tap on the rows on the display to see more details and dose history.

To view dose history stored in the dosimeter, connect the dosimeter to a computer running Dose Viewer (**C** in Figure “2.” on page 2). See the *Dose Viewer manual*, found in the help menu in Dose Viewer.

If you require further information about product training, contact your local RaySafe representative.

CAUTION! *The Real-time Dosimeter measures scattered radiation and should neither be positioned in the primary X-ray field nor be used for diagnostic measurements.*

WARNING! *RaySafe i3 does not replace any X-ray protection equipment.*

WARNING! *RaySafe i3 is not sterile.*

WARNING! *RaySafe i3 is not a replacement for legal dosimeters.*

Maintenance

The Real-time Dosimeter has a replaceable battery. When a red light is blinking on the side of the dosimeter, the battery needs to be replaced. A service technician can do this, using the instructions in the service manual.

The Real-time Display has no serviceable parts.

For cleaning of dosimeter and display, use a damp cloth and mild detergent.

WARNING! *Do not modify the Real-time Display.*

Specifications

Physical

	Real-time Dosimeter	Real-time Display
Weight	34 g (1.2 oz)	1.2 kg (42 oz)
Dimensions	40 × 58 × 17 mm (1.6 × 2.3 × 0.7 in)	300 × 250 × 60 mm (11.8 × 9.8 × 2.4 in)
Power	Replaceable battery: CR2450 Li/MnO ₂	External supply: Astec DPS53-M, GTM96600-6012-T3 or FW7556M/12

Radiology

Operational quantity	H _p (10)
Energy dependence	< 25 % (N-series, 40 – 150 kV)
Temperature dependence	< 5 % (18 – 26 °C) < 25 % (15 – 18 °C, 26 – 35 °C)
Response time	< 1 s (> 100 µSv/h), < 5 s (< 100 µSv/h)
Integration interval	1 s
Detection limit	< 30 µSv/h
Dose rate uncertainty (continuous radiation)	10 % or 10 µSv/h (40 µSv/h – 150 mSv/h), 20 % (150 – 300 mSv/h), 40 % (300 – 500 mSv/h)
The measured dose rate is monotonically increasing up to 1 Sv/h.	
Dose reproducibility	10 % or 1 µSv
Maximum lifetime dose	10 Sv
Typical angular response (relative to 0°)	Vertical: <i>Figure “3.” on page 2</i> Horizontal: <i>Figure “4.” on page 2</i>
Typical energy response	<i>Figure “5.” on page 2</i>
Real-time Display output unit	Sv or rem

Radio

	Real-time Dosimeter	Real-time Display
Frequency	868.3 MHz, 918.3 MHz or 927.9 MHz, depending on purchased configuration (see label on product)	
Transmitter power output	3.2 mW	1 mW
Communication range	Typically 5 – 15 m, depending on the local environment	

Climatic conditions

Temperature	+15 – +35 °C (for operation) –25 – +70 °C (for storage)
Temperature rate of change	< 0.5 °C/minute
Relative humidity	20 – 80 % (for operation) 5 – 95 % (for storage)
IP code	Real-time Dosimeter: IP20 Real-time Display: IP21

For safe usage

Always power off the dosimeter before air transport, to make sure that the radio communication is completely turned off.

If there is any risk of liquid spilling on the Real-time Display, turn it off by disconnecting the power cord.

Security and privacy

Dosimeter serial number as well as name and settings entered via Dose Viewer are sent unencrypted via radio to the Real-time Display.

Symbols on labels

Symbol	Meaning
	Manufacturer
	Date of manufacture
	Article number
	Serial number
	Consult the instructions for use.
	Product complies with essential requirements applicable to European directives
	Electrical safety compliance marking, together with listing number (“Appendix” on page 15)
	Radio certification number for US market (“Appendix” on page 15)
	Radio certification number for Canadian market (“Appendix” on page 15)
	Radio certification symbol and number for Brazilian market (“Appendix” on page 15)
	Conforms to relevant Australian Safety and EMC standards.
	Radio certification symbol and number for Japanese market
	In the European Union the WEEE directive and label indicates that when discarding this product, it must be sent to separate collection facilities for recovery and recycling. For advice and information, contact RaySafe.
	None of the China RoHS restricted substances are present above permitted levels.

Learn more

- Real-time Dosimeter quick guide under the lid of the dosimeter packaging
- Dose Viewer manual via the Help button in Dose Viewer
- Installation and Service manual delivered with the RaySafe i3 system

Service and support contacts

For information about service and support, visit www.raysafe.com.

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APPENDIX

MET

MET Laboratories, Inc. Certification covers UL 60950-1/CSA C22.2 No. 60950-1. MET Laboratories, Inc. has not evaluated this device for reliability or efficacy of its intended functions.

FCC

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This portable device with its antenna is excluded from SAR testing since the output power is below the SAR test exclusion threshold. The antenna used for this device must not be co-located or operating in conjunction with any other antenna or transmitter.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

IC/ISED

This Device complies with Industry Canada License-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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1256-14-5254

Sustainability

For declaration of regulated substances,
see <https://www.raysafe.com/about-us/sustainability>

EU DECLARATION OF CONFORMITY

Name / address manufacturer: Unfors RaySafe AB
Uggledalsv. 29
SE-427 40 Billdal
Sweden

This declaration of conformity is issued under the sole responsibility of the manufacturer.

Product name: RaySafe i3 868.3 MHz

The product consists of the
following sub-products:

RaySafe Real-time Display 868.3 MHz	1606036
RaySafe Real-time Dosimeter 868.3 MHz	1601027
RaySafe Real-time Hub 868.3 MHz	1601034
RaySafe Video Unit	4760036

Product Options/Accessories: Accessories as described in the
accompanying labeling documentation

The object of the declaration described above is in conformity with the following Directives, specified per product in table below:

• <i>Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment. Amended by Directive (EU) 2015/863 of the European Parliament and of the Council of 31 March 2015 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.</i>	RaySafe Real-time Display RaySafe Real-time Dosimeter RaySafe Real-time Hub RaySafe Video Unit
• <i>Directive 2014/30/EU on the harmonisation of the laws of the Member States relating to electromagnetic compatibility</i>	RaySafe Video Unit
• <i>Directive 2014/53/EU of the European Parliament and of the Council of 16 April 2014 on the harmonisation of the laws of the Member States relating to the making available on the market of radio equipment and repealing Directive 1999/5/EC</i>	RaySafe Real-time Display RaySafe Real-time Dosimeter RaySafe Real-time Hub

The product(s) listed above have been tested in a typical configuration as described in the Manufacturer's accompanying documentation and are fully compliant with the standards listed below. Additionally, the products listed above have been designed, manufactured, tested, and found to be compatible with the devices and accessories described by the manufacturer in the devices accompanying documentation.

Date of Issue: 2019-12-09

Signature:
(signed for and on
behalf of RaySafe)

Name: Daniel Carlsson
Title: Manager RA/QA

Place of issue: Billdal, Sweden

The products as mentioned above to which this EU Declaration of conformity is applicable, has been developed, designed and found to be in conformity pursuant to the European harmonized standards specified hereafter

Aspect	Harmonized Standard	Title	Affected product
Safety (RED Directive)	EN60950-1:2006 /A11:2009/A1:2010/ A12:2011/A2:2013	Information technology equipment - Safety – Part 1: General requirements	RaySafe Real-time Display RaySafe Real-time Dosimeter RaySafe Real-time Hub (RaySafe Video Unit is not a radio product but it's still approved against this standard)
Radio (RED Directive)	ETSI EN 300 220-1 V3.1.1	Short Range Devices (SRD) operating in the frequency range 25 MHz to 1 000 MHz; Part 1: Technical characteristics and methods of measurement	RaySafe Real-time Display RaySafe Real-time Dosimeter RaySafe Real-time Hub
Radio (RED Directive)	ETSI EN 300 220-2 V3.1.1	Short Range Devices (SRD) operating in the frequency range 25 MHz to 1 000 MHz; Part 2: Harmonised Standard covering the essential requirements of article 3.2 of Directive 2014/53/EU for nonspecific radio equipment	RaySafe Real-time Display RaySafe Real-time Dosimeter RaySafe Real-time Hub
EMC (RED and EMC Directive)	ETSI EN 301 489-1 V2.1.1	Electromagnetic Compatibility (EMC)= standard for radio equipment and services; Part 1: Common technical requirements; Harmonised Standard covering the essential requirement of article 3.1(b) of Directive 2014/53/EU and the essential requirements of article 6 of Directive 2014/30/EU	RaySafe Real-time Display RaySafe Real-time Dosimeter RaySafe Real-time Hub RaySafe Video Unit
EMC (RED and EMC Directive)	ETSI EN 301 489-3 V2.1.1	Electromagnetic Compatibility (EMC) standard for radio equipment and services; Part 3: Specific conditions for Short-Range Devices (SRD) operating on frequencies between 9 kHz and 246 GHz; Harmonised Standard covering the essential requirements of article 3.1(b) of Directive 2014/53/EU	RaySafe Real-time Display RaySafe Real-time Dosimeter RaySafe Real-time Hub RaySafe Video Unit
Environmental (RoHS Directive)	EN 50581:2012	Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances	Real-time Display Real-time Dosimeter Real-time Hub RaySafe Video Unit